



PRODUCT SPECIFICATION

BEARD SNOOD WHITE NON-STERILE

FaceGuard  5

PRODUCT INFORMATION

PRODUCT DESCRIPTION

White latex free beard snood with elasticated ear loops

PRODUCT FEATURES

One Size with elasticated ear loops
Manufactured from spun-bonded polypropylene

PACKAGING

Supplied in a polythene bag inside a cardboard box

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems

SHELF LIFE

3 years

STORAGE CONDITIONS

Store in cool, well ventilated area

RECYCLING INFORMATION

Component Name:	Component Material:	Disposal:	Guidance:
Beard Snood	Polypropylene Fabric with elasticated ear loops	Recyclable	Remove elastic ear lops
Polythene bag	LDPE	Recyclable	Check With Local Authority, Remove Any Labels
Packaging Tape	PE Tape	Non-recyclable	Remove From Cardboard Box
Box	Cardboard	Recyclable	Remove Tape
Box Label	Coated Paper With a Hot Glue Adhesive	Non-recyclable Individually	Leave On Cardboard Box
Product Label	Coated Paper With a Hot Glue Adhesive	Non-recyclable	Remove From Polybag

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty
MNWH01ZZZ	Beard Snood White Non-sterile	One Size	100	500