



## PRODUCT SPECIFICATION

### SINGLE USE CLEANROOM COMPATIBLE SOCK

TuffGuard 

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Fully compliant with the requirements of EU GMP Annex 1

A 100% synthetic fibre, single use cleanroom sock

The sock is 'sock shaped' for increased wearer comfort

##### PRODUCT FEATURES

White 98% Polyester / 2% Elastaine Sock

Sock shaped

One Size

##### PARTICULATE TESTING

The polyester socks have been tested for particulate levels using the Helmke Drum method described in the standard IEST RP CC 003.4

##### Applicable standards

The socks meet Category II of the standard

##### PACKAGING

100 Socks (50 pairs) are packed in a single sealed polythene bag

##### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems

##### STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions

## RECYCLING INFORMATION

Component Name	Component Material	Disposal	Guidance
Sock	98% Polyester / 2% Elastaine	Non-recyclable	Mixed fibre so check local recycling facilities
Box	Cardboard	Recyclable	Check with local authority, remove any labels and tape
Polybag	LDPE	Recyclable	Check with local authority, remove any labels and tape

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
SNWH50001	Single Use Cleanroom Compatible Sock	One Size	100 (50 pairs)	100 (50 pairs)	400 x 310 x 110 mm	2.00 KG