



## PRODUCT SPECIFICATION

### WHITE OVERSHOE WITH TEXTURED PVC SOLE NON-STERILE



#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Overshoe with white textured pvc sole

##### PRODUCT FEATURES

White  
Stitched seams with sole  
Weight of shoe cover - 50g+380g

##### PACKAGING

Primary packaged in polythene packaging

##### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems

##### STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions

##### MATERIAL AND PERFORMANCE

Material - Microporous polypropylene with polyethylene laminate fabric + Polyvinylchloride (PVC) soles  
Overshoe size  
Back Length - 12.5cm/tolerance (1cm)  
Width - 42cm/tolerance (1cm)  
Sole Length - 33.5cm/tolerance (1cm)

#### ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty
SNWH87016		White Overshoe with textured PVC sole Non-sterile	One Size	20	20 packs of 20